

SAM[®]

JUNCTIONAL TOURNIQUET

TO CONTROL DIFFICULT BLEEDS IN THE INGUINAL & AXILLA AREAS

TO IMMOBILIZE A PELVIC FRACTURE

The SAM[®] Junctional Tourniquet for hemorrhage control is designed to control bleeding in areas where standard tourniquets would not be effective, such as with IED/Blast injuries or high level amputations.

- Multiple Indications
 - Inguinal hemorrhage
 - Axilla hemorrhage
 - Pelvic immobilization
- Lightweight (1 lb., 1 oz.) (488g)
- Easy to use (4 step application)
- < 25 second application time
- FDA 510(k) cleared and CE marked

A patented buckle maintains correct base force to eliminate slack. “Click” provides clear feedback to confirm when to secure strap to Velcro for correct application.

Large, repositionable and/or removable Target Compression Devices (TCD) ensure easy targeting. A relief valve prevents over-inflation.

Easy to apply. Quick Application (typically under 25 seconds with training).



Inflates up to 3.2”.



Clinically proven to provide the prescribed safe and effective force for pelvic stabilization.

Removable hand pump quickly inflates the Target Compression Devices (TCD) to achieve hemorrhage control.

